

## GENERAL FLUOROSCOPIC EQUIPMENT DATA FORMS

MEDICAL/DENTAL X-RAY EQUIPMENT DATA																																																	
NAVMED 6470/4 (7-80)		REPORT SYMBOL MED 6470-15																																															
<b>1. FACILITY IDENTIFICATION</b>																																																	
a. FACILITY NAME	b. UIC																																																
c. MAILING ADDRESS	d. BUILDING	e. ROOM																																															
<b>2. STATUS OF THE EQUIPMENT</b> (INDICATE IF EQUIPMENT IS IN USE OR THE REASON FOR NOT BEING IN USE). <div style="display: flex; justify-content: space-between;"> <div> <input type="checkbox"/> IN USE  <input type="checkbox"/> NOT IN USE         </div> <div> <input type="checkbox"/> TO BE REPAIRED  <input type="checkbox"/> CANNOT BE REPAIRED         </div> <div> <input type="checkbox"/> STORED IN GOOD WORKING CONDITION  <input type="checkbox"/> OTHER         </div> </div>																																																	
<b>3. X-RAY EQUIPMENT IDENTIFICATION</b> <div style="display: flex; justify-content: space-between; margin-top: 10px;"> <div style="width: 40%;">           a. PLANT ACCOUNT NUMBER            b. YEAR EQUIPMENT WAS MANUFACTURED            c. INSTALLATION DATE OF EQUIPMENT            d. X-RAY EQUIPMENT IS CERTIFIED: YES ____ NO ____         </div> <div style="width: 55%;"> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr><td style="height: 20px;"></td></tr> <tr><td style="height: 20px;"></td></tr> <tr><td style="height: 20px;"></td></tr> </table> </div> </div> <div style="margin-top: 20px;"> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 30%;">e. COMPONENT</th> <th style="width: 30%;">f. MANUFACTURER</th> <th style="width: 20%;">g. MODEL</th> <th style="width: 20%;">h. SERIAL NUMBER</th> </tr> </thead> <tbody> <tr><td>1) CONTROL CONSOLE</td><td></td><td></td><td></td></tr> <tr><td>2) X-RAY TABLE</td><td></td><td></td><td></td></tr> <tr><td>3) X-RAY TUBE ASSEMBLY</td><td></td><td></td><td></td></tr> <tr><td style="padding-left: 20px;">TUBE #1 HOUSING</td><td></td><td></td><td></td></tr> <tr><td style="padding-left: 20px;">TUBE #1 INSERT</td><td></td><td></td><td></td></tr> <tr><td style="padding-left: 20px;">TUBE #1 COLLIMATOR</td><td></td><td></td><td></td></tr> <tr><td> </td><td></td><td></td><td></td></tr> <tr><td>IMAGE INTENSIFIER</td><td></td><td></td><td></td></tr> <tr><td></td><td></td><td></td><td></td></tr> <tr><td></td><td></td><td></td><td></td></tr> </tbody> </table> </div> <div style="text-align: right; margin-top: 10px;"> <input type="checkbox"/> CONTINUED ON SEPARTE SHEET         </div>						e. COMPONENT	f. MANUFACTURER	g. MODEL	h. SERIAL NUMBER	1) CONTROL CONSOLE				2) X-RAY TABLE				3) X-RAY TUBE ASSEMBLY				TUBE #1 HOUSING				TUBE #1 INSERT				TUBE #1 COLLIMATOR								IMAGE INTENSIFIER											
e. COMPONENT	f. MANUFACTURER	g. MODEL	h. SERIAL NUMBER																																														
1) CONTROL CONSOLE																																																	
2) X-RAY TABLE																																																	
3) X-RAY TUBE ASSEMBLY																																																	
TUBE #1 HOUSING																																																	
TUBE #1 INSERT																																																	
TUBE #1 COLLIMATOR																																																	
IMAGE INTENSIFIER																																																	
<b>4. TYPE OF X-RAY EQUIPMENT (CHECK AS MANY AS APPROPRIATE)</b> <div style="display: flex; justify-content: space-between;"> <div> <input type="checkbox"/> RADIOGRAPHIC  <input type="checkbox"/> FLUOROSCOPIC  <input type="checkbox"/> COMBINATION R/F         </div> <div> <input type="checkbox"/> FIXED  <input type="checkbox"/> MOBILE  <input type="checkbox"/> OTHER _____         </div> <div> <input type="checkbox"/> DENTAL INTRAORAL  <input type="checkbox"/> DENTAL PANOGRAPHIC         </div> </div>																																																	
<b>5. GENERATOR (CHECK ONE)</b> <div style="display: flex; justify-content: space-between;"> <div> <input type="checkbox"/> AUTORECTIFIED  <input type="checkbox"/> SINGLE PHASE HALF WAVE  <input type="checkbox"/> SINGLE PHASE FULL WAVE         </div> <div> <input type="checkbox"/> THREE PHASE  <input type="checkbox"/> CAPACITOR DISCHARGE  <input type="checkbox"/> OTHER (SPECIFY) _____         </div> <div>           MAXIMUM mA ____ mA             MAXIMUM kVp ____ kVp         </div> </div>																																																	
<b>6. ASSOCIATED EQUIPMENT (CHECK AS MANY AS APPROPRIATE)</b> <div style="display: flex; justify-content: space-between;"> <div> <input type="checkbox"/> AUTOMATIC EXPOSURE CONTROL SYSTEM (LIKE PHOTOTIMER)  <input type="checkbox"/> SPOT FILM DEVICE         </div> <div> <input type="checkbox"/> OTHER _____         </div> <div> <input type="checkbox"/> PHOTOSPOT CAMERA  <input type="checkbox"/> IMAGE INTENSIFIER         </div> </div>																																																	
<b>7. USE (CHECK ONE)</b> <div style="display: flex; justify-content: space-between;"> <div> <input type="checkbox"/> GENERAL RADIOGRAPHY  <input type="checkbox"/> CHEST RADIOGRAPHY  <input type="checkbox"/> HEAD RADIOGRAPHY         </div> <div> <input type="checkbox"/> MAMMOGRAPHY  <input type="checkbox"/> TOMOGRAPHY  <input type="checkbox"/> UROLOGY STUDIES         </div> <div> <input type="checkbox"/> OTHER (SPECIFY) _____         </div> </div>																																																	
<b>8. DATE OF LAST RADIATION PROTECTION SURVEY</b>  DATE: <input type="checkbox"/> UNKNOWN	<b>9. THIS EQUIPMENT REPLACED EQUIPMENT WITH PLANT ACCOUNT NUMBER.</b>  <input type="checkbox"/> UNKNOWN																																																
10. REPORTED BY:	REVIEWED BY:	DATE:																																															
TITLE:																																																	

GENERAL REQUIREMENTS AND PERFORMANCE TESTS FOR FLUOROSCOPIC EQUIPMENT			
NAVMED 6470/7 (10-99)		REPORT SYMBOL MED 6470-10	
<b>I. FACILITY IDENTIFICATION</b>			
1. FACILITY NAME	2. UIC	3. BUILDING	4. ROOM
5. MAILING ADDRESS			
<b>II. X-RAY EQUIPMENT IDENTIFICATION</b>			
1. X-RAY MACHINE: CERTIFIED YES ____ NO ____			
2. X-RAY TUBE:	a. MANUFACTURER:		
	b. MODEL:		
	c. SERIAL NUMBER:		
3. IMAGE DEVICE:	a. TYPE:		
	b. SERIAL NUMBER:		
<b>III. OPERATIONAL AND RADIATION SAFETY CHARACTERISTICS</b>			
1. PROTECTION DEVICES		YES	NO
a. LEAD GLOVES AND APRONS: 0.25 mm LEAD EQUIVALENCE OR MORE.			
b. TOWER APRON; 0.25 mm LEAD EQUIVALENCE OR MORE.			
c. BUCKY SLOT COVER; 0.25 mm LEAD EQUIVALENCE OR MORE.			
d. EXPOSURE SWITCH REQUIRES CONTINUOUS PRESSURE TO OPERATE.			
e. FLUOROSCOPY TUBE DOES NOT PRODUCE X-RAYS UNLESS THE IMAGE DEVICE ASSEMBLY IS IN POSITION INTERCEPTING THE X-RAY BEAM.			
f. THE X-RAY FIELD DOES NOT EXTEND BEYOND THE IMAGE DEVICE ASSEMBLY.			
g. SHUTTERS CLOSE OR REDUCE THE X-RAY FIELD TO NO MORE THAT 5 x 5 cm. AT MAXIMUM SID.			
h. CUMULATIVE "ON TIME" (5 MIN MAX) TIMER IS PROVIDED AND IS MANUALLY RESET.			
i. AUDIBLE SIGNAL INDICATES COMPLETION OF PRESET CUMULATIVE TIME.			
j. HIGH LEVEL CONTROL, IF AVAILABLE, REQUIRES THE OPERATOR TO APPLY CONTINUOUS PRESSURE AND AN AUDIBLE SIGNAL INDICATES CONTROL IS BEING USED.			
k. CERTIFIED EQUIPMENT: LENGTH OR WIDTH OF X-RAY FIELD IN THE PLANE OF IMAGE RECEPTOR DOES NOT EXCEED VISIBLE AREA BY MORE THAN 3% OF SID.			
l. SPOT FILM DEVICE IS FUNCTIONING SATISFACTORLY.			
2. REMARKS-			

GENERAL REQUIREMENTS AND PERFORMANCE TESTS FOR FLUOROSCOPIC EQUIPMENT (CONT)

NAVMED 6470/7 (10-99)

REPORT SYMBOL MED 6470-10

IV. PERFORMANCE TESTS

1. KILOVOTAGE ACCURACY

a. kVp SETTING		b. kVp DETERMINED		c. ACCURACY	
1)	50 kVp		kVp		
2)	60 kVp		kVp		
3)	70 kVp		kVp		
4)	80 kVp		kVp		
5)	90 kVp		kVp		
6)	100 kVp		kVp		
7)	110 kVp		kVp		
8)	120 kVp		kVp		
9)			kVp		
10)			kVp		

d. kV CHECKING DEVICE USED:

1) TYPE:

2) MODEL:

3) SERIAL NUMBER:

2. TYPICAL ENTRANCE EXPOSURE RATE

IMAGE SIZE:	Phantom:	b. kVP	c. mA	d. R/MIN
1. MANUAL: 2. AUTOMATIC BRIGHTNESS CONTROL: 3. HIGH LEVEL CONTROL 4. PULSE FLUOROSCOPY 5. CINE				
1. MANUAL: 2. AUTOMATIC BRIGHTNESS CONTROL: 3. HIGH LEVEL CONTROL 4. PULSE FLUOROSCOPY 5. CINE				
1. MANUAL: 2. AUTOMATIC BRIGHTNESS CONTROL: 3. HIGH LEVEL CONTROL 4. PULSE FLUOROSCOPY 5. CINE				
1. MANUAL: 2. AUTOMATIC BRIGHTNESS CONTROL: 3. HIGH LEVEL CONTROL 4. PULSE FLUOROSCOPY 5. CINE				

SURVEYOR:

DATE:

GENERAL REQUIREMENTS AND PERFORMANCE TESTS FOR FLUOROSCOPIC EQUIPMENT (CON'T)																									
NAVMED 6470/7 (10-99)			REPORT SYMBOL MED 6470-10																						
<b>IV. PERFORMANCE TESTS (CONTINUED)</b>																									
3. MAXIMUM ENTRANCE EXPOSURE RATE																									
IMAGE SIZE:	Phantom:		b. kVP	c. mA	d. R/MIN																				
1. MANUAL:																									
2. AUTOMATIC BRIGHTNESS CONTROL:																									
3. HIGH LEVEL CONTROL																									
4. PULSE FLUOROSCOPY																									
5. CINE																									
MAXIMUM ENTRANCE EXPOSURE RATE SHOULD NOT BE MORE THAN:																									
10 R/MIN - UNCERTIFIED EQUIPMENT																									
5 R/MIN - CERTIFIED EQUIPMENT AND ALL EQUIPMENT WITH HIGH LEVEL CONTROL, WHEN IT IS NOT IN USE.																									
4. RADIATION TRANSMITTED THROUGH PRIMARY BARRIER																									
a. Technique used to obtain maximum entrance exposure rate			b. MEASUREMENT (mR/hr)																						
MODE	kVp	mA																							
INSTRUMENT USED																									
a. TYPE	b. MODEL	c. SERIAL NUMBER	d. CALIBRATION DATE																						
RADCAL MDH			DATE:																						
5. BEAM QUALITY																									
a. MEASUREMENTS		b. I.I. SIZE	c. Mode	d. kVP	e. mA																				
					f. OTHER																				
<div style="display: flex; justify-content: space-around;"> <div style="text-align: center;"> <p>ADDED FILTRATION</p> <table border="1" style="margin: auto;"> <tr><td>1</td><td></td></tr> <tr><td>2</td><td></td></tr> <tr><td>3</td><td></td></tr> <tr><td>4</td><td></td></tr> <tr><td>5</td><td></td></tr> </table> </div> <div style="text-align: center;"> <p>EXPOSURE</p> <table border="1" style="margin: auto;"> <tr><td>1</td><td></td></tr> <tr><td>2</td><td></td></tr> <tr><td>3</td><td></td></tr> <tr><td>4</td><td></td></tr> <tr><td>5</td><td></td></tr> </table> </div> </div>		1		2		3		4		5		1		2		3		4		5		<p>g. HVL: _____ mm Al. EQUIVALENT.</p> <p>h. EQUIPMENT COMPLIES THE HVL REQUIREMENTS:</p> <p style="text-align: center;">YES    ____    NO    ____</p>			
1																									
2																									
3																									
4																									
5																									
1																									
2																									
3																									
4																									
5																									
6. MINIMUM SSD																									
MINIMUM SOURCE TO SKIN DISTANCE _____ CM.			MODE OF DETERMINATION:		TRIANGULATION OR TAPE MEASURE																				
7. BEAM LIMITATION																									
a. MINIMUM FIELD SIZE _____			MAX SID:	LARGEST I.I. SIZE:																					
b. MAXIMUM FIELD SIZE _____																									
8. FLUOROSCOPIC DISPLAY SIZE																									
	VERTICAL AXIS (cm)	HORIZONTAL AXIS (cm)	SID																						
VIEWING FIELD (TV)																									
RADIATION FIELD (FILM)																									
ALIGNMENT DIFF (% OF SID)																									
SUM OF ALIGNMENT DIFFERENCES EXPRESSED AS A PERCENTAGE OF THE SID			<div style="border: 1px solid black; width: 80px; height: 20px; margin: auto;"></div>																						
9. BEAM CENTRAL ALIGNMENT																									
MINIMUM SID	MAXIMUM I.I. SIZE																								
IS BEAM AXIS/I.I. MISALIGNMENT LESS THAN OR EQUAL TO 1.5 DEGREES			YES	NO																					
			<div style="border: 1px solid black; width: 80px; height: 20px; margin: auto;"></div>	<div style="border: 1px solid black; width: 80px; height: 20px; margin: auto;"></div>																					

GENERAL REQUIREMENTS AND PERFORMANCE TESTS FOR FLUOROSCOPIC EQUIPMENT (CONT)									
NAVMED 6470/7 (10-99)						REPORT SYMBOL MED 6470-10			
<b>IV. PERFORMANCE TESTS (CONTINUED)</b>									
10. SPATIAL DISTORTION									
FLUOROSCOPIC MODE: MANUAL OR ABC					LARGEST IMAGE INTENSIFIER SIZE				
a. PIN-CUSHION DISTORTION			YES		NO				
b. "S-ING" DISTORTION			YES		NO				
11. HIGH CONTRAST RESOLUTION									
					NUMBER OF LINE PAIRS VISIBLE AT EDGE/CENTER			MAX SID	
FLUOROSCOPIC MODE (I.I. SIZE)									
MANUAL MODE									
AUTOMATIC BRIGHTNESS CONTROL									
HIGH LEVEL CONTROL									
PULSED FLUOROSCOPY									
CINE									
SPOT FILMING (MECHANICAL)									
SPOT FILMING (DIGITAL)									
12. LOW CONTRAST SENSITIVITY									
					HOLE SIZE VISIBLE			MAX SID	
FLUOROSCOPIC MODE (I.I. SIZE)									
MANUAL MODE									
AUTOMATIC BRIGHTNESS CONTROL									
HIGH LEVEL CONTROL									
PULSED FLUOROSCOPY									
CINE									
SPOT FILMING (MECHANICAL)									
SPOT FILMING (DIGITAL)									
IS 3.1 mm TEST HOLE AT 2% CONTRAST VISIBLE FOR ALL MODES EXCEPT PULSE FLUORO						YES / NO			
14. AUTOMATIC EXPOSURE CONTROL (AEC) SYSTEM									
SID:		FILMSCREEN COMBINATION:			FILM SIZE:				
a. OPTICAL DENSITY EVALUATION:									
kVp	mA	Detector Cell	AEC Setting	Phantom Thickness	Image Number	Elapsed mAs	Measured OD		
b. OUTPUT REPRODUCIBILITY									
Imaging Mode		kVp	mA	AEC Setting	SID				
Detector Cell		Output (mR)		Reading 1	Reading 2	Reading 3	Mean	-5%	+5%
							#DIV/0!	#DIV/0!	#DIV/0!
ARE ALL VALUES WITHIN +/- 5% OF MEAN				YES / NO					
c. BACK-UP TIMER									
Imaging Mode		kVp	mA	Lead Thickness		AEC Setting	SID		
Detector Cell		Elapsed mAs		DOES THE BEAM TERMINATE PRIOR TO 600 mAs					
				YES			NO		
d. kVp COMPENSATION									
mA		Phantom Thickness		AEC Setting					
kVp		Elapsed mAs		Optical Density					
						Density Range			
						Do all densities lie within			
						+/- 3% of baseline density			
						determined in 13a. Yes / No			

REPORT SYMBOL MED 6470-10

	IMAGE GENERATED BY IMAGE HOLD	DIGITAL SPOT FILM	
STEP 1 OD			SID OF 30 cm
STEP 2 OD			MAX I.I. SIZE
STEP 3 OD			
STEP 4 OD			Plot density as a function of
STEP 5 OD			wedge thickness for both
STEP 6 OD			images and compare

## GENERAL REQUIREMENTS AND PERFORMANCE TESTS FOR FLUOROSCOPIC EQUIPMENT (CONT)

NAVMED 6470/7 (10-99)

REPORT SYMBOL MED 6470-10

## IV. PERFORMANCE TESTS (CONTINUED)

## CONCLUSIONS

	YES	NO
a. MAXIMUM ENTRANCE EXPOSURE RATE DOES NOT EXCEED 10 RMIN UNDER NORMAL NORMAL CONDITIONS (NOT EQUIPPED WITH HIGH LEVEL CONTROL).		
b. MAXIMUM ENTRANCE EXPOSURE RATE DOES NOT EXCEED 5 RMIN WHEN EQUIPPED WITH HIGH LEVEL CONTROL.		
c. THE SOURCE TO SKIN DISTANCE FOR UNDER TABLE FLUOROSCOPIC TUBE IS NOT LESS THAN 12".		
d. RADIATION TRANSMITTED THROUGH PRIMARY BARRIER DOES NOT EXCEED 1 mR/HR AT 10 cm FROM PRIMARY BARRIER REAR SURFACE.		
e. THE HALF-VALUE LAYER WAS DETERMINED TO BE _____ IN mm Al EQUIVALENTS.		
f. EQUIPMENT COMPLIES WITH HALF-VALUE REQUIREMENTS.		
g. kVp DETERMINED TO BE WITHIN +/- 5 % OF NOMINAL SETTING OR READOUT.		
h. PERFORMANCE TESTS ON IMAGE INTENSIFIER WERE SATISFACTORY.		
i. AEC OUTPUT REPRODUCIBILITY FOR SPOT FILM DEVICE IS WITHIN +/- 5% OF MEAN.		
j. OPTICAL DENSITY RANGES FOR SPOT FILM AEC EVALUATION TESTS WITHIN ACCEPTABLE LIMITS.		

SURVEYOR:

DATE: